

TECO DIAGNOSTICS

1268 N. Lakeview Ave.
Anaheim, CA 92807
1-800-222-9880

CARBON DIOXIDE (CO₂) REAGENT SET

INTENDED USE

For manual or automated quantitative determination of Carbon Dioxide content in human plasma or serum at 340 nm.

INTRODUCTION

Carbon Dioxide (CO₂) in serum or plasma exists primarily as dissolved CO₂ and bicarbonate anion (HCO₃⁻).¹ The plasma CO₂ content is decreased in metabolic acidosis and respiratory alkalosis, whereas the level is increased in metabolic alkalosis and respiratory acidosis.² In pathologic conditions such as in diabetes mellitus, glomerulonephritis, pyloric obstruction, diarrhea, etc., acidosis or alkalosis could be anticipated.³ Therefore, determination of plasma CO₂ content as part of an electrolyte profile can help establish to a degree, the anticipated change in the patient.

Plasma CO₂ content can be measured manometrically or using an ion selective electrode or by spectrophotometric procedures.⁴ This CO₂ reagent measures CO₂ content enzymatically and the procedure is a modification of the method of Forrester et al.⁵

PRINCIPLES OF THE PROCEDURE



Phosphoenol pyruvate carboxylase (PEPC) catalyzes the reaction between phosphoenol pyruvate and carbon dioxide (bicarbonate) to form oxalacetate and phosphate ion. Oxalacetate is reduced to malate with simultaneous oxidation of an equimolar amount of reduced nicotinamide adenine dinucleotide (NADH) to NAD; the reaction is catalyzed by malate dehydrogenase (MDH). This results in a decrease in absorbance at 340 nm that is directly proportional to CO₂ concentration in the sample.

REAGENT COMPOSITION

1. When reconstituted according to the direction, CO₂ reagent contains: PEP 1.8 mM, Magnesium Sulfate 10 mM, NADH 0.40 mM, MDH (porcine) 1,250 U/L, PEPC (microbial) 200 U/L, Sodium Oxamate 2.5 mM, Buffer, (pH 7.0), Non-reactive fillers and stabilizers with Sodium Azide 0.1% as preservative.
2. CO₂ Standard contains 30 mmol/L of Sodium bicarbonate in an aqueous solution.

WARNINGS AND PRECAUTIONS

1. CO₂ reagent is for "In Vitro Diagnostic Use."
2. Normal precautions exercised in handling laboratory reagents should be followed.
3. CO₂ reagent contains sodium azide that is toxic if ingested and which may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide accumulation.
4. Use only CO₂-free water for reconstitution.
5. Do not pipette by mouth to avoid CO₂ contamination from the expired air.

REAGENT PREPARATION

Reconstitute CO₂ Reagent with the volume of CO₂-free water, indicated on the vial label. Mix by inversion 5-6 times. Do not shake. CO₂ standard is ready to use.

Note: Avoid contamination of reagents with CO₂. Do not blow into pipette, since breath contains a high content of CO₂. Do not let bottles open unnecessarily, since CO₂ from air can contaminate reagent. Keep container tightly stoppered.

STORAGE AND STABILITY

1. Store the unopened vial of CO₂ Reagent and CO₂ standard refrigerated (2 - 8°C). Reagent and standard are stable until the expiration date shown on the labels.
2. Reconstituted reagents stored in the closed tightly stoppered bottle are stable for 1 day at room temperature (18 - 26°C) and 7 days refrigerated (2 - 8°C).

REAGENT DETERIORATION

Do not use the reagent if the absorbance of CO₂ reagent when measure against water in a 1-cm light path is less than 1.000 at 340 nm, or dry reagent exhibit caking due to possible moisture penetration or do not dissolve completely.

SPECIMEN COLLECTION AND STORAGE

1. It is recommended that specimen collection be carried out in accordance with NCCLS document M29-T
2. Serum or heparinized plasma may be assayed. Ideally, venous blood should be collected and handled anaerobically.
3. Plasma and serum, after prompt separation from cells or clot, should be kept tightly stoppered.
4. CO₂ content of blood is stable for 1 hour when stored at 2-4°C under anaerobic conditions.^{3,6} Alkalinized serum can stand, without loss of CO₂, for as long as 4 hours in open cups.⁷

INTERFERING SUBSTANCES

The major interference in this assay is CO₂ from air or from the breath of the analyst. Certain drugs and other substances are also known to influence blood CO₂ levels. See Young DS.⁸

MATERIALS PROVIDED

Carbon Dioxide Reagent and Carbon Dioxide Standard.

MATERIALS REQUIRED, BUT NOT PROVIDED

Spectrophotometer, capable of accurately measuring absorbances at 340 nm; Cuvettes with optical properties suitable for use at 340 nm; Pipetting devices for the accurate delivery of volumes required for this assay; Timer; Incubator or water bath, capable of maintaining 37°C.

PROCEDURE

1. Prepare CO₂ Reagent according to REAGENT PREPARATION.
2. Label tubes "Blank", "Standard", "Controls", "Patients", etc.
3. Pipette 1.0 mL Carbon Dioxide Reagent into each tube.
4. Incubate all tubes for 3 minutes at 37°C.
5. Set spectrophotometer wavelength at 340 nm, temperature to 37°C.

- Pipette 5 µL (0.005 mL) of water, standard, and sample to the cuvette labeled "Blank", "Standard" and "Patients", respectively.
- Mix gently by inversion and incubate for 5 minutes.
- Read and record absorbance (Abs.) of all cuvettes at 340 nm.

PROCEDURE NOTES

If the samples is lipemic, icteric, or hemolyzed, a sample blank must be prepared by adding 5 µl of sample to 1.0 mL of saline. The absorbance of this blank is subtracted from the absorbance of the respective test and the corrected absorbance is then used in the calculations.

QUALITY CONTROL

The integrity of the reaction should be monitored by use of normal and abnormal controls with known CO₂ values for each assay.

CALIBRATION

CO₂ Standard provided is standardized for using in manual procedure calculation. The concentration of CO₂ should be validated by comparison with commercially available standard. For automated instrument, refer to appropriate application specification.

CALCULATIONS

Determine CO₂ content of sample as follows:

$$\text{CO}_2 \text{ Content of Sample (mmol/L)} = \frac{\text{Abs. blank} - \text{Abs. sample}}{\text{Abs. blank} - \text{Abs. standard}} \times \text{concentration of standard}$$

Example:

$$\begin{aligned} \text{Abs. sample} &= 0.95 \\ \text{Abs. standard} &= 0.82 \\ \text{Abs. blank} &= 1.66 \\ \text{Concentration of Standard} &= 30 \text{ mmol/L} \end{aligned}$$

$$\begin{aligned} \text{CO}_2 \text{ (mmol/L)} &= \frac{1.66 - 0.95}{1.66 - 0.82} \times 30 \\ &= 25.4 \text{ mmol/L} \end{aligned}$$

PROCEDURAL LIMITATIONS

CO₂ reagent is linear to 40 mmol/L. If CO₂ content is greater than 40 mmol/L, dilute sample with an equal volume of CO₂-free isotonic saline and reassay. Multiply the result by 2 to compensate for dilution. Lipemic, icteric, or hemolyzed sample will influence the result. Sample blank must be prepared. Please refer to PROCEDURE NOTES section for the preparation of sample blank.

EXPECTED VALUES

23-34 mmol/L.⁶

It is strongly recommended that each laboratory establish an expected range characteristic for the local population.

PERFORMANCE

- Linearity:** 0-40 mmol/L.
- Sensitivity:** Based on an instrument resolution of Abs.= 0.001, this procedure has a sensitivity of 0.03 mmol/L.
- Comparison:**

Manual procedure: 40 serum samples with CO₂ range from 7.0-30.4 mmol/L were assayed. The study against similar method yielded the correlation coefficient of 0.987 and linear regression equation $y = 1.02x - 0.42$

Automated procedure: 70 serum samples with CO₂ range from 0.7-29.7 mmol/L were assayed. The study against similar method yielded the correlation coefficient of 0.99 and linear regression equation $y = 1.03x - 0.69$

- Precision:** The following studies were performed according to NCCLS guideline #EP5-T. Each control was subjected to 20 replicated assays for each study.

Within Run	Control	Mean (mmol/l)	S.D. (mmol/l)	C.V. (%)
Manual procedure	I	34.0	0.6	1.8
	II	16	0.7	4.1
Automated procedure	III	23.13	0.50	2.2
	IV	15.67	0.47	3.0


Run to Run	Control	Mean (mmol/l)	S.D. (mmol/l)	C.V. (%)
Manual procedure	I	33.4	1.4	4.3
	II	15.7	0.7	4.4
Automated procedure	III	23.29	2.00	8.6
	IV	14.79	1.60	10.8

REFERENCES

- Clinical Chemistry, LA Kaplan, AJ Pesce, Editors, CV Mosby Company, St. Louis (MO), p. 1056 (1984).
- Hydrogen Ion Concentration in Body Fluids. IN Contarow and Trumper Clinical Biochemistry, 7th ed., AL Latner, Editor, Saunders, Philadelphia, p.399 (1975).
- Tietz NW, Pruden EL, Siggaard-Andersen O: Electrolytes, Blood Gases and Acid-Base Balance. IN Textbook of Clinical Chemistry, NW Tietz, Editor, Saunders, Philadelphia, p. 1188 (1986).
- Natelson S: Routine use of ultramicro methods in the clinical laboratory. Am J Clin Pathol 21:1153, (1951).
- Forrester RL, Wataji LJ, Silverman DA, Pierre KJ; Enzymatic method for the determination of CO₂ in serum. Clin Chem 22: 243, (1976).
- Ibbott FA, LaGanga TS, Gin JB, Inkpen JA: Blood pH, CO₂ and O₂. IN Clinical Chemistry- Principles and Technics, RJ Henry, DC Cannon, JW Winkelman, Editors, Harper & Row, Hagerstown (MD), p 755 (1974).
- Gambino SR, Schreiber H; The measurement of CO₂ content with the Auto Analyzer: A comparison with 3 standard methods and a description of a new method (alkalinization) for preventing loss of CO₂ from open cups. Am J Clin Pathol 45:406, (1966).
- Young DS: Effects of Drugs on clinical Laboratory Tests, 3rd ed., AACCC Press, (1990).
- Wolf PL: Interpretation of Biochemical Multitest Profiles. An Analysis of 100 Important Conditions. Masson Publishing, New York, p 33,56 (1977).

C525: 02/2018

Manufactured by:



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands